

# Fusion Vision System with Epi Fluorescence

View opaque samples via reflected light.

## Features:

- Populate the fluorescence Block with Single or Multiband filter sets
- Variable wavelengths of light input via LED light source
- Choose from 7:1 and 12.5:1 zoom ranges
- Built in focusing and fine focus adjustment
- Specify use with DSLR or C-mount cameras
- Variety of mounting and clamping options
- Motorised or Manual operation
- Incorporate LWD objectives
- Also use via a microscope with our 38mm ISO coupler system

